

AKAS DVT GARMENT

Helps to prevent Deep Vein Thrombosis

- Tailored for Indian skins.
- Breathable fabric to prevent sweat accumulation
- 100% made in India and supporting local manufacturing
- Available in a range of size to provide perfect fit



Feature	HELLO 111	HELLO 123
Compression Mode	Uniform Only	Sequential & Uniform
Weight	2 kg	2.5 kg
Battery Backup	Up to 35 hours	Up to 24 hours
Pressure Range	Thigh & Calf: 30-60 mmHg	Thigh & Calf: 30-60 mmHg
	Foot: 120-150 mmHg	Foot: 120-150 mmHg
Charging Time	5 hrs (normal), 3 hrs (fast)	5 hrs (normal), 3 hrs (fast)

Thigh Garment
 L - 71.1 to 91.4 cm
 M - 55.9 to 71.7 cm
 S - 40.6 to 55.9 cm

Calf Garment
 L - 53 to 66 cm
 M - 35.6 to 53 cm
 S - Up to 35.6 cm

Foot Garment
 M - US 9, UK 8.5

Manufactured by:



SP 127, 1st Main Road, Ambattur Industrial Estate, Chennai- 600058
 M : +91-98403 79116 | E : sales@akasinfusions.com | W : www.akasinfusions.com

Hello 111

Bye Bye DVT



Say **Hello** to Healing, and **Goodbye** to DVT

Up to 60% of DVT cases can be prevented. HELLO aids in safe healing through smart compression technology, which enhances circulation and quietly and effectively lowers the risk of DVT.

Hello 111

Our DVT pump promotes circulation and helps prevent DVT, lymphedema, and chronic venous insufficiency in at-risk or immobile patients.



✦ Responsive Therapy

Automatically adjusts pressure and cycle timing based on patient body type and garment for safe, targeted therapy.

✦ Therapy Made Easy

One-touch startup ensures therapy begins quickly and easily, reducing stress for both patients and caregivers



24 Hours Battery Backup

✦ Backup Beyond Limits

Delivers 24 hours of continuous therapy—ensuring uninterrupted DVT protection throughout recovery without the need to recharge.

✦ Fast Charging

The first ever DVT pump with fast charging—just 1 hour of charge delivers up to 6 hours of therapy.



Robust Build



Long Lasting Battery



Bed Hook Support



All Angle Alarm



Bigger Display



2kg weight



Standby mode